

TNNI1 Antibody, FITC conjugated



PACO50064

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Gene ID:

TNNI1

Uniprot

P19237

Synonyms:

Troponin I, slow skeletal muscle (Troponin I, slow-twitch isoform), TNNI1

Immunogen:

Recombinant Human Troponin I, slow skeletal muscle protein (2-173AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A