

XPO1 Antibody, Biotin conjugated



PACO50101

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

XPO1

Uniprot

O14980

Synonyms:

Exportin-1 (Exp1) (Chromosome region maintenance 1 protein homolog), XPO1, CRM1

Immunogen:

Recombinant Human Exportin-1 protein (625-845AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A