

CIDEB Antibody



PACO50154

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Activates apoptosis.

Gene ID:

CIDEB

Uniprot

Q9UHD4

Synonyms:

Cell death activator CIDE-B (Cell death-inducing DFFA-like effector B), CIDEB


Immunogen:

Recombinant Human Cell death activator CIDE-B protein (80-145AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO50154 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung tissue using PACO50154 at dilution of 1:100.