

CIDEB Antibody, HRP conjugated



PACO50155

Product Information

Size:	Protein Background:
50ug	Activates apoptosis.
Reactivity:	Gene ID:
Human	CIDEB
Source:	Uniprot
Rabbit	Q9UHD4
Isotype:	Synonyms:
IgG	Cell death activator CIDE-B (Cell death-inducing DFFA-like effector B), CIDEB
Applications:	Immunogen:
ELISA	Recombinant Human Cell death activator CIDE-B protein (80-145AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A