

PRKCZ Antibody



PACO50182

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Gene ID:

PRKCZ

Uniprot

Q05513

Synonyms:

Protein kinase C zeta type (EC 2.7.11.13) (nPKC-zeta), PRKCZ, PKC2


Immunogen:

Recombinant Human Protein kinase C zeta type protein (494-575AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO50182 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).