

Allergen Fel d 4 Antibody



PACO50198

Product Information

Size:

50ug

Reactivity:

Cat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be a pheromone carrier. Acts as a kairomone, detected by the prey vomeronasal organ and inducing fear reactions in mice.

Gene ID:**Uniprot**

Q5VFH6

Synonyms:

Allergen Fel d 4 (allergen Fel d 4)

Immunogen:

Recombinant Cat Allergen Fel d 4 protein (16-186AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A