

Major allergen Bla g 1.0101 Antibody



PACO50202

Product Information

Size:

50ug

Reactivity:

Blattella germanica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:**Gene ID:**

Major allergen Bla g 1.0101

Uniprot

Q9UAM5


Synonyms:**Immunogen:**

Recombinant Blattella germanica Major allergen Bla g 1.0101 protein (1-412AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: Major allergen Bla g 1.0101 antibody at 2.8 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 62 kDa. Observed band size: 62 kDa.