

# Allergen Bla g 4 Antibody



PACO50206

---

## Product Information

**Size:**

50ug

**Reactivity:**

Blattella germanica

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable ligand-binding protein.

**Gene ID:**

**Uniprot**

P54962

**Synonyms:**

Allergen Bla g 4 (Allergen Bla g IV) (allergen Bla g 4) (Fragment)

**Immunogen:**

Recombinant Blattella germanica Allergen Bla g 4 protein (13-182AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A