

# Major pollen allergen Aln g 1 Antibody



PACO50210

---

## Product Information

**Size:**

50ug

**Reactivity:**

Alnus glutinosa

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

**Uniprot**

P38948

**Synonyms:**

Major pollen allergen Aln g 1 (Allergen Aln g I) (allergen Aln g 1)

**Immunogen:**


Recombinant Alnus glutinosa Major pollen allergen Aln g 1 protein (2-160AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: Major pollen allergen Aln g 1 antibody at 3 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 34 kDa. observed band size: 34 kDa.