

**Major pollen allergen Aln g 1 Antibody, HRP
conjugated
PACO50211**



Product Information

Size:

50ug

Reactivity:

Alnus glutinosa

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P38948

Synonyms:

Major pollen allergen Aln g 1 (Allergen Aln g I) (allergen Aln g 1)

Immunogen:

Recombinant Alnus glutinosa Major pollen allergen Aln g 1 protein (2-160AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A