

ALK Antibody, Biotin conjugated



PACO50237

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ALK

Uniprot

Q9UM73

Synonyms:

ALK tyrosine kinase receptor (EC 2.7.10.1) (Anaplastic lymphoma kinase) (CD antigen CD246), ALK

Immunogen:

Recombinant Human ALK tyrosine kinase receptor protein (1402-1604AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A