

PACO50274

---

## Product Information

**Size:**

50ug

**Reactivity:**

Klebsiella pneumoniae

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

mrkA

**Uniprot**

P12267

**Synonyms:**

Fimbrial subunit type 3, mrkA

**Immunogen:**

Recombinant Klebsiella pneumoniae Fimbrial subunit type 3 protein (23-202AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A