

bla Antibody



PACO50282

Product Information

Size:

50ug

Reactivity:

Pseudomonas aeruginosa

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Hydrolyzes both carbenicillin and oxacillin.

Gene ID:

bla

Uniprot

P14489

Synonyms:

Beta-lactamase OXA-10 (EC 3.5.2.6) (Beta-lactamase PSE-2), bla, oxa10 pse2

Immunogen:

Recombinant Pseudomonas aeruginosa Beta-lactamase OXA-10 protein (21-266AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A