

p37 Antibody



PACO50294

Product Information

Size:

50ug

Reactivity:

Mycoplasma hyorhinitis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

P37 is part of a high-affinity transport system.

Gene ID:

p37

Uniprot

P15363

Synonyms:

High affinity transport system protein p37, p37

Immunogen:

Recombinant Mycoplasma hyorhinitis High affinity transport system protein p37 protein (24-403AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A