

PACO50298

Product Information

Size:

50ug

Reactivity:

Chlamydia pneumoniae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Prevents misfolding and promotes the refolding and proper assembly of unfolded polypeptides generated under stress conditions.

Gene ID:

groL

Uniprot

P31681

Synonyms:

60 kDa chaperonin (GroEL protein) (Protein Cpn60), groL, groEL mopA

Immunogen:

Recombinant Chlamydia pneumoniae 60 kDa chaperonin protein (1-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A