

Product Information

Size:

50ug

Reactivity:

Puumala virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Encapsidates the genome protecting it from nucleases. The encapsidated genomic RNA is termed the nucleocapsid (NC) and serves as template for transcription and replication.

Gene ID:

N

Uniprot

P41268

Synonyms:

Nucleoprotein (Nucleocapsid protein) (Protein N), N

Immunogen:

Recombinant Puumala virus Nucleoprotein protein (1-433AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A