

BETVII Antibody



PACO50318

Product Information

Size:

50ug

Reactivity:

Betula pendula

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds to actin and affects the structure of the cytoskeleton. At high concentrations, profilin prevents the polymerization of actin, whereas it enhances it at low concentrations. By binding to PIP2, it inhibits the formation of IP3 and DG.

Gene ID:

BETVII

Uniprot

P25816

Synonyms:

Profilin-1 (Allergen Bet v II) (Pollen allergen Bet v 2) (allergen Bet v 2), BETVII

Immunogen:

Recombinant Betula pendula Profilin-1 protein (2-133AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A