

Beta-lytic metalloendopeptidase Antibody, HRP conjugated
PACO50359



Product Information

Size:

50ug

Reactivity:

Lysobacter enzymogenes

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P00801

Synonyms:

Beta-lytic metalloendopeptidase (EC 3.4.24.32) (Beta-lytic protease)

Immunogen:

Recombinant Lysobacter enzymogenes Beta-lytic metalloendopeptidase protein (1-178AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A