

PACO50422

---

## Product Information

**Size:**

50ug

**Reactivity:**

Penaeus monodon

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:2000-1:10000

**Protein Background:**

arginine kinase activity, phosphorylation.

**Gene ID:**

AK

**Uniprot**

C7E3T4

**Synonyms:**

Arginine kinase (AK) (EC 2.7.3.3) (allergen Pen m 2), AK

**Immunogen:**


Recombinant Penaeus monodon Arginine kinase protein (2-356AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: AK antibody at 2.8 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 55 kDa. Observed band size: 55 kDa.