
Product Information

Size:

50ug

Reactivity:

Neosartorya fumigata

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

intracellular, IgE binding.

Gene ID:

aspf3

Uniprot

O43099

Synonyms:

Peroxiredoxin Asp f3 (Prx) (EC 1.11.1.15) (Thioredoxin peroxidase) (TPx) (allergen Asp f 3), aspf3


Immunogen:

Recombinant Neosartorya fumigata Peroxiredoxin Asp f3 protein (1-168AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: aspf3 antibody at 3.25 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 35 kDa. Observed band size: 35 kDa.