

PON2 Antibody, Biotin conjugated



PACO50489

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Capable of hydrolyzing lactones and a number of aromatic carboxylic acid, esters. Has antioxidant activity. Is not associated with high density lipoprotein. Prevents LDL lipid peroxidation, reverses the oxidation of mildly oxidized LDL, and inhibits the ability of MM-LDL to induce monocyte chemotaxis.

Gene ID:

PON2

Uniprot

Q15165

Synonyms:

Serum paraoxonase/arylesterase 2 (PON 2) (EC 3.1.1.2) (EC 3.1.1.81) (Aromatic esterase 2) (A-esterase 2) (Serum arylalkylphosphatase 2), PON2

Immunogen:

Recombinant Human Serum paraoxonase/arylesterase 2 protein (32-97AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A