

# CAMK1 Antibody, Biotin conjugated



PACO50505

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CAMK1

**Uniprot**

Q14012

**Synonyms:**

Calcium/calmodulin-dependent protein kinase type 1 (EC 2.7.11.17) (CaM kinase I) (CaM-KI) (CaM kinase I alpha) (CaMKI-alpha), CAMK1

**Immunogen:**

Recombinant Human Calcium/calmodulin-dependent protein kinase type 1 protein (166-265AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A