

MYBPH Antibody, HRP conjugated



PACO50535

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds to myosin; probably involved in interaction with thick myofilaments in the A-band.

Gene ID:

MYBPH

Uniprot

Q13203

Synonyms:

Myosin-binding protein H (MyBP-H) (H-protein), MYBPH

Immunogen:

Recombinant Human Myosin-binding protein H protein (3-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A