

DGKZ Antibody, FITC conjugated



PACO50540

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Displays a strong preference for 1,2-diacylglycerols over 1,3-diacylglycerols, but lacks substrate specificity among molecular species of long chain diacylglycerols. Isoform 2 but not isoform 1 regulates RASGRP1 activity.

Gene ID:

DGKZ

Uniprot

Q13574

Synonyms:

Diacylglycerol kinase zeta (DAG kinase zeta) (EC 2.7.1.107) (Diglyceride kinase zeta) (DGK-zeta), DGKZ, DAGK6

Immunogen:

Recombinant Human Diacylglycerol kinase zeta protein (11-152AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A