

**Guanine nucleotide-binding protein subunit beta-like protein Antibody, FITC conjugated**  
**PACO50612**



---

**Product Information**

**Size:**

50ug

**Reactivity:**

Glycine max

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

**Uniprot**

Q39836

**Synonyms:**

Guanine nucleotide-binding protein subunit beta-like protein

**Immunogen:**

Recombinant Glycine max Guanine nucleotide-binding protein subunit beta-like protein (1-325AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A