

Pollen allergen Phl p 5b Antibody



PACO50622

Product Information

Size:

50ug

Reactivity:

Phleum pRatense

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Has ribonuclease activity. May be involved in host-pathogen interactions.

Gene ID:**Uniprot**

Q40963

Synonyms:

Pollen allergen Phl p 5b (Allergen Phl p Vb) (allergen Phl p 5b) (Fragment)


Immunogen:

Recombinant Phleum pratense Pollen allergen Phl p 5b protein (20-284AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: Pollen allergen Phl p 5b antibody at 2.7µg/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 43 kDa. observed band size: 43 kDa.