

PACO50642

Product Information

Size:

50ug

Reactivity:

Paracoccus denitrificans

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Nitrous-oxide reductase is part of a bacterial respiratory system which is activated under anaerobic conditions in the presence of nitrate or nitrous oxide.

Gene ID:

nosZ

Uniprot

Q51705

Synonyms:

Nitrous-oxide reductase (EC 1.7.2.4) (N(2)OR) (N2O reductase), nosZ

Immunogen:

Recombinant Paracoccus denitrificans Nitrous-oxide reductase protein (58-553AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A