

hlgC Antibody



PACO50658

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Toxin that seems to act by forming pores in the membrane of the cell. Has a hemolytic and a leucotoxic activity.

Gene ID:

hlgC

Uniprot

Q6GE13

Synonyms:

Gamma-hemolysin component C, hlgC

Immunogen:

Recombinant Staphylococcus aureus Gamma-hemolysin component C protein (30-315AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A