

Major pollen allergen Art v 1 Antibody



PACO50742

Product Information

Size:

50ug

Protein Background:

Reactivity:

Artemisia vulgaris

Gene ID:

Source:

Rabbit

Uniprot

Q84ZX5

Isotype:

IgG

Synonyms:

Major pollen allergen Art v 1 (Defensin-like protein 1) (allergen Art v 1)

Applications:

ELISA, WB

Immunogen:

Recombinant Artemisia vulgaris Major pollen allergen Art v 1 protein (25-132AA).


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: Major pollen allergen Art v 1 antibody at 3.4 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 27 kDa. observed band size: 27 kDa.