

TPH2 Antibody, Biotin conjugated



PACO50761

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, indolalkylamine biosynthetic process.

Gene ID:

TPH2

Uniprot

Q8IWU9

Synonyms:

Tryptophan 5-hydroxylase 2 (EC 1.14.16.4) (Neuronal tryptophan hydroxylase) (Tryptophan 5-monoxygenase 2), TPH2, NTPH

Immunogen:

Recombinant Human Tryptophan 5-hydroxylase 2 protein (44-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A