

SULF2 Antibody, HRP conjugated



PACO50763

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Exhibits arylsulfatase activity and highly specific endoglucosamine-6-sulfatase activity. It can remove sulfate from the C-6 position of glucosamine within specific subregions of intact heparin.

Gene ID:

SULF2

Uniprot

Q8IWU5

Synonyms:

Extracellular sulfatase Sulf-2 (hSulf-2) (EC 3.1.6), SULF2, KIAA1247

Immunogen:

Recombinant Human Extracellular sulfatase Sulf-2 protein (463-660AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A