

ABCA5 Antibody



PACO50862

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May play a role in the processing of autolysosomes.

Gene ID:

ABCA5

Uniprot

Q8WWZ7

Synonyms:

ATP-binding cassette sub-family A member 5, ABCA5, KIAA1888


Immunogen:

Recombinant Human ATP-binding cassette sub-family A member 5 protein (1230-1355AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO50862 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO50862 at dilution of 1:100.