

NEUROD6 Antibody, HRP conjugated



PACO50871

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Activates E box-dependent transcription in collaboration with TCF3/E47. May be a trans-acting factor involved in the development and maintenance of the mammalian nervous system. Transactivates the promoter of its own gene.

Gene ID:

NEUROD6

Uniprot

Q96NK8

Synonyms:

Neurogenic differentiation factor 6 (NeuroD6) (Class A basic helix-loop-helix protein 2) (bHLHa2) (Protein atonal homolog 2), NEUROD6, ATOH2 BHLHA2

Immunogen:

Recombinant Human Neurogenic differentiation factor 6 protein (170-314AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A