

VPS33A Antibody



PACO50890

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

VPS33A

Uniprot

Q96AX1

Synonyms:

Vacuolar protein sorting-associated protein 33A (hVPS33A), VPS33A


Immunogen:

Recombinant Human Vacuolar protein sorting-associated protein 33A protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO50890 at dilution of 1:100.



Immunofluorescent analysis of U251 cells using PACO50890 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).