

# P2RX4 Antibody, Biotin conjugated



PACO50973

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Receptor for ATP that acts as a ligand-gated ion channel. This receptor is insensitive to the antagonists PPADS and suramin.

**Gene ID:**

P2RX4

**Uniprot**

Q99571

**Synonyms:**

P2X purinoceptor 4 (P2X4) (ATP receptor) (Purinergic receptor), P2RX4

**Immunogen:**

Recombinant Human P2X purinoceptor 4 protein (91-176AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A