

YIPF3 Antibody



PACO51010

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200, IF:1:50-1:200

Protein Background:

Involved in the maintenance of the Golgi structure. May play a role in hematopoiesis.

Gene ID:

YIPF3

Uniprot

Q9GZM5

Synonyms:

Protein YIPF3 (Killer lineage protein 1) (Natural killer cell-specific antigen KLIP1) (YIP1 family member 3) [Cleaved into: Protein YIPF3, 36 kDa form III], YIPF3, C6orf109 KLIP1


Immunogen:

Recombinant Human Protein YIPF3 protein (2-148AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4


Product Images



Western Blot. Positive WB detected in: Mouse heart tissue, Mouse liver tissue, Mouse brain tissue, Mouse stomach tissue. All lanes: YIPF3 antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 39 kDa. Observed band size: 39 kDa.



Immunofluorescent analysis of PC-3 cells using PACO51010 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human colon cancer using PACO51010 at dilution of 1:100.