

X Antibody



PACO51026

Product Information

Size:

50ug

Reactivity:

Hepatitis B virus genotype D subtype ayw

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

X

Uniprot

Q9QMI3

Synonyms:

Protein X (HBx) (Peptide X) (pX), X

Immunogen:

Recombinant Hepatitis B virus genotype D subtype ayw Protein X protein (1-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A