

TRIM39 Antibody, FITC conjugated



PACO51100

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

E3 ubiquitin-protein ligase. May facilitate apoptosis by inhibiting APC/C-Cdh1-mediated poly-ubiquitination and subsequent proteasome-mediated degradation of the pro-apoptotic protein MOAP1.

Gene ID:

TRIM39

Uniprot

Q9HCM9

Synonyms:

E3 ubiquitin-protein ligase TRIM39 (EC 2.3.2.27) (RING finger protein 23) (RING-type E3 ubiquitin transferase TRIM39) (Testis-abundant finger protein) (Tripartite motif-containing protein 39), TRIM39, RNF23 TFP

Immunogen:

Recombinant Human E3 ubiquitin-protein ligase TRIM39 protein (310-416AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A