

# KLF2 Antibody, Biotin conjugated



PACO51229

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcription factor that binds to the CACCC box in the promoter of target genes such as HBB/beta globin or NOV and activates their transcription.

**Gene ID:**

KLF2

**Uniprot**

Q9Y5W3

**Synonyms:**

Krueppel-like factor 2 (Lung krueppel-like factor), KLF2, LKLF

**Immunogen:**

Recombinant Human Krueppel-like factor 2 protein (1-355AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A