

SEC23IP Antibody, HRP conjugated



PACO51243

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in the organization of endoplasmic reticulum exit sites. Specifically binds to phosphatidylinositol 3-phosphate (PI(3)P), phosphatidylinositol 4-phosphate (PI(4)P) and phosphatidylinositol 5-phosphate (PI(5)P).

Gene ID:

SEC23IP

Uniprot

Q9Y6Y8

Synonyms:

SEC23-interacting protein (p125), SEC23IP

Immunogen:

Recombinant Human SEC23-interacting protein (779-989AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A