

ANXA13 Antibody



PACO51314

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Involved in vesicular traffic to the apical plasma membrane (isoform B).

Gene ID:

ANXA13

Uniprot

Q29471

Synonyms:

Annexin A13 (Annexin XIII) (Annexin-13) (Intestine-specific annexin) (ISA), ANXA13, ANX13


Immunogen:

Recombinant Dog Annexin A13 protein (2-316AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: ANXA13 antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 36 kDa. observed band size: 36 kDa.