

# ARG1 Antibody, FITC conjugated



PACO51324

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, cytosol, extracellular exosome, nucleus, arginase activity, arginine catabolic process, urea cycle.

**Gene ID:**

ARG1

**Uniprot**

P05089

**Synonyms:**

Arginase-1 (EC 3.5.3.1) (Liver-type arginase) (Type I arginase), ARG1

**Immunogen:**

Recombinant Human Arginase-1 protein (1-322AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A