

ASNS Antibody, HRP conjugated



PACO51331

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, asparagine synthase (glutamine-hydrolyzing) activity, protein homodimerization activity, asparagine biosynthetic process, cellular amino acid, biosynthetic process, cellular response to glucose starvation, glutamine metabolic process, negative regulation of apoptotic process, PERK-mediated unfolded protein response, positive regulation of mitotic cell cycle.

Gene ID:

ASNS

Uniprot

P08243

Synonyms:

Asparagine synthetase [glutamine-hydrolyzing] (EC 6.3.5.4) (Cell cycle control protein TS11) (Glutamine-dependent asparagine synthetase), ASNS, TS11

Immunogen:

Recombinant Human Asparagine synthetase [glutamine-hydrolyzing] protein (364-551AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A