

# AVPR2 Antibody, HRP conjugated



PACO51343

---

## Product Information

**Size:**

50ug

**Reactivity:**

Pig

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Receptor for arginine vasopressin. The activity of this receptor is mediated by G proteins which activate adenylate cyclase. Involved in renal water reabsorption.

**Gene ID:**

AVPR2

**Uniprot**

P32307

**Synonyms:**

Vasopressin V2 receptor (V2R) (AVPR V2) (Antidiuretic hormone receptor) (Renal-type arginine vasopressin receptor), AVPR2

**Immunogen:**

Recombinant Pig Vasopressin V2 receptor protein (292-370AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A