

PACO51370

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Reversible hydration of carbon dioxide. Its role in saliva is unknown.

Gene ID:

CA6

Uniprot

P23280

Synonyms:

Carbonic anhydrase 6 (EC 4.2.1.1) (Carbonate dehydratase VI) (Carbonic anhydrase VI) (CA-VI) (Salivary carbonic anhydrase) (Secreted carbonic anhydrase), CA6


Immunogen:

Recombinant Human Carbonic anhydrase 6 protein (17-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human salivary gland tissue using PACO51370 at dilution of 1:100.