

EIF2AK1 Antibody, Biotin conjugated



PACO51481

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EIF2AK1

Uniprot

Q9BQI3

Synonyms:

Eukaryotic translation initiation factor 2-alpha kinase 1 (EC 2.7.11.1) (Heme-controlled repressor) (HCR) (Heme-regulated eukaryotic initiation factor eIF-2-alpha kinase) (Heme-regulated inhibitor) (Hemin-sensitive initiation factor 2-alpha kinase), EIF2AK1, HRI KIAA1369

Immunogen:

Recombinant Human Eukaryotic translation initiation factor 2-α kinase 1 protein (1-86AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A