

PACO51502

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays an important role in the regulation of embryonic development, cell proliferation, and cell differentiation.

Gene ID:

fgf3

Uniprot

P48802

Synonyms:

Fibroblast growth factor 3 (FGF-3) (Heparin-binding growth factor 3) (HBGF-3), fgf3, fgf-3

Immunogen:

Recombinant Zebrafish Fibroblast growth factor 3 protein (19-256AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A