

FUCA1 Antibody, FITC conjugated



PACO51520

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Alpha-L-fucosidase is responsible for hydrolyzing the alpha-1,6-linked fucose joined to the reducing-end N-acetylglucosamine of the carbohydrate moieties of glycoproteins.

Gene ID:

FUCA1

Uniprot

P04066

Synonyms:

Tissue alpha-L-fucosidase (EC 3.2.1.51) (Alpha-L-fucosidase I) (Alpha-L-fucoside fucosylase 1) (Alpha-L-fucosidase 1), FUCA1

Immunogen:

Recombinant Human Tissue alpha-L-fucosidase protein (32-158AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A