

Gfap Antibody, HRP conjugated



PACO51539

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

GFAP, a class-III intermediate filament, is a cell-specific marker that, during the development of the central nervous system, distinguishes astrocytes from other glial cells.

Gene ID:

Gfap

Uniprot

P47819

Synonyms:

Glial fibrillary acid, c protein (GFAP), Gfap

Immunogen:

Recombinant Rat Glial fibrillary acid, c protein (1-430AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A