

GZMK Antibody, FITC conjugated



PACO51600

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

serine-type peptidase activity.

Gene ID:

GZMK

Uniprot

P49863

Synonyms:

Granzyme K (EC 3.4.21) (Fragmentin-3) (Granzyme-3) (NK-tryptase-2) (NK-Tryp-2), GZMK, TRYP2

Immunogen:

Recombinant Human Granzyme K protein (27-264AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A